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Feeding chutes and funnels

To help particles flow appropriately during the test process, the BeVision D2 prepares several customization options for feeding chutes and the feeding funnel.

The customizable coating material for the funnel and chute can effectively prevent particles from sticking to the wall. The funnel size and the chute size are also customizable based on requirements, permitting the BeVision D2 to hold an appropriate volume of samples.

Adjustable feeding guides

The feeding guide is used to orient the particle flow in the focusing range of the objective lens, thus the quality of captured images and accuracy of measurements are ensured.

Based on the approximate particle size, the width of feeding guide openings can be adjusted to achieve a balance between the feeding rate and the quality of images.

Objective lens

Depending on the approximate particle size and the width of the particle size distribution, the BeVision D2 offers 5 optional objective lenses to ensure measurement accuracy.

Each optical lens has its own efficient working range, to accurately measure particle size and shape of a great variety of samples. According to ISO 13322, it is recommended to select a working range covering the whole size distribution range of the sample.

Objective lens	1 μm	10 μm	100 μm	1,000 μm	10,000 μm
2x lens	3.75	11.25	39.75	500	
1x lens	7.5	22.5	67.5	1800	
0.5x lens	15	45	135	3600	
0.3x lens	25	75	225	6000	
0.132x lens	56.81	170.43	511.29	13655	

■ able to detect particles, rough size measurement
 ■ size and rough shape measurement
 ■ size and shape measurement in ISO-area

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